

Date: August 07, 2026

Author: Software Engineering

System: Apex

Subsystem: Microsoft Win Pro Embedded 10 LTSB 2016

Subject: Microsoft Win Pro Embedded 10 LTSB 2016 Patch release July 2026

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSB 2016 systems.

Note: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the "Sentinel Run-Time 7.80.0.0" which can be located in the Windows Control Panel "Program and Features". If needed you can install the patch by ordering "APEX Service Dongle Patch Release" (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2014:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for July 2026:

KB5099542
KB4589210

KB5099535

KB5090407

List of all patches/KB's validated on the Apex System

KB5090407	KB5061010	KB5036899
KB5099535	KB5060954	KB5035962
KB5099542	KB5058524	KB5035855
KB4589210	KB5058383	KB5034862
KB5094122	KB5055661	KB5034767
KB5094141	KB5055521	KB5034127
KB5087537	KB5054833	KB5034119
KB5088064	KB890830	KB5033904
KB5082198	KB5054006	KB5033688
KB5082089	KB5053594	KB5033373
KB5068404	KB5052006	KB5033371
KB5078938	KB5050109	KB5032391
KB5075902	KB5049993	KB5032197
KB5075999	KB5048671	KB5032196
KB5073447	KB5046612	KB5031984
KB5073722	KB5044293	KB5031362
KB5071543	KB5043124	KB5031361
KB5065222	KB5043051	KB5030504
KB5068864	KB5042749	KB5030214
KB5070247	KB5041773	KB5030213
KB5066584	KB5041576	KB5029931
KB5066836	KB5040562	KB5029247
KB5065427	KB5040434	KB5029242
KB5065687	KB5039334	KB5029185
KB5063871	KB5039214	KB5028169
KB5062560	KB5037763	KB5028168
KB5062799	KB5037016	KB5027222

KB5027219	KB5016058	KB5005112
KB5027131	KB5015811	KB5005043
KB5026363	KB5015808	KB5005030
KB5026362	KB5014702	KB5004244
KB5025229	KB5014692	KB5004238
KB5025228	KB5014553	KB5003711
KB5023788	KB5014026	KB5003646
KB5023702	KB5013952	KB5003638
KB5023697	KB5013941	KB5003243
KB5022840	KB5013641	KB5003197
KB5022838	KB5012647	KB5003171
KB5022511	KB5012596	KB5002542
KB5022289	KB5012170	KB5002520
KB5022286	KB5012128	KB5002497
KB5021237	KB5011570	KB5002464
KB5021235	KB5011503	KB5002406
KB5021126	KB5011495	KB5001407
KB5021045	KB5010359	KB5001404
KB5020866	KB5010351	KB5001402
KB5020627	KB5009557	KB5001347
KB5019966	KB5009546	KB5001342
KB5019964	KB5008873	KB5001078
KB5018419	KB5008218	KB5000859
KB5018411	KB5008207	KB5000822
KB5017396	KB5007206	KB5000803
KB5017315	KB5007192	KB5000685
KB5017305	KB5006672	KB4601393
KB5017095	KB5006669	KB4601392
KB5016884	KB5005698	KB4601345
KB5016623	KB5005573	KB4601318
KB5016622	KB5005568	KB4601060

KB4598480	KB4558997	KB4514366
KB4598243	KB4557397	KB4512937
KB4598230	KB4556813	KB4512578
KB4593226	KB4552924	KB4512577
KB4592440	KB4551853	KB4512574
KB4589210	KB4550994	KB4512517
KB4589208	KB4550929	KB4511553
KB4587735	KB4549949	KB4509091
KB4586830	KB4549947	KB4507460
KB4586793	KB4540723	KB4506998
KB4583465	KB4540670	KB4505225
KB4580346	KB4539571	KB4503308
KB4580325	KB4538461	KB4498947
KB4577668	KB4537764	KB4494440
KB4577667	KB4537759	KB4494175
KB4577586	KB4535680	KB4494174
KB4577015	KB4534273	KB4493470
KB4576750	KB4534271	KB4489882
KB4571694	KB4532691	KB4487026
KB4570333	KB4530715	KB4485447
KB4570332	KB4530689	KB4346087
KB4566424	KB4525236	KB4132216
KB4565912	KB4523205	KB4091664
KB4565511	KB4523204	KB4022619
KB4565349	KB4521862	KB4018073
KB4562562	KB4521858	KB2565063
KB4562561	KB4520724	KB2467173
KB4561616	KB4519998	
KB4561608	KB4519338	
KB4561600	KB4516115	
KB4558998	KB4516044	