

Date: September 04, 2026

Author: Software Engineering

System: Apex

Subsystem: Win Pro Embedded 10 LTSC 2019

Subject: Microsoft Win Pro Embedded 10 LTSC 2019 Patch release August 2026

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSC 2019 systems.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2019:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated:

<u>Date</u>	<u>List of KB's</u>		
July 2026	KB5100989	KB5099538	KB5090407
August 2026	KB5120238	KB5120698	

List of all patches/KB's validated on the Apex System before July 2026:

KB5100989	KB5044277	KB5023702
KB5099538	KB5044016	KB5022840
KB5090407	KB5043050	KB5022511
KB5094123	KB5042749	KB5022286
KB5087061	KB5041578	KB5021237
KB5087538	KB5040430	KB5021045
KB5082123	KB5039879	KB5020866
KB5082413	KB5039217	KB5019964
KB5068404	KB5037765	KB5018419
KB5078752	KB5036896	KB5017396
KB5075904	KB5036604	KB5017315
KB5073723	KB5035849	KB5016623
KB5071544	KB5034768	KB5015811
KB5065222	KB5034127	KB5014692
KB5068791	KB5033904	KB5014164
KB5066143	KB5033688	KB5013941
KB5066586	KB5033371	KB5013641
KB5065428	KB5032196	KB5012647
KB5063877	KB5031984	KB5012170
KB5062557	KB5031361	KB5012128
KB5060531	KB5030214	KB5011503
KB5058392	KB5029931	KB5010351
KB5055519	KB5029247	KB5009557
KB5054833	KB5028168	KB5008873
KB890830	KB5027222	KB5008218
KB5053596	KB5027131	KB5007206
KB5052000	KB5026362	KB5006672
KB5050008	KB5025229	KB5005698
KB5049608	KB5005568	KB5003711
KB5048661	KB5005112	KB5002542
KB5046615	KB5005030	KB5002520
KB5002497	KB5002223	KB4577586
KB5002464	KB4589210	KB4535680
KB5002406	KB4589208	KB4494174
KB5002369	KB4583465	KB4132216
KB5002323	KB4583462	KB4022619
KB4018073		

