



Date: September 9, 2026

Author: Software Engineering

System: BREVERA / CORLUMINA

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release
August 2026

Purpose:

To release the list of validated Microsoft critical patches on the BREVERA / CORLUMINA System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated:

<u>Date</u>	<u>KB's validated:</u>		
July 2026	KB5100989	KB5099538	KB5090354
August 2026	KB5120698	KB5120238	

List of all patches/KB's validated on the BREVERA / CORLUMINA System before July 2026:

KB5100989	KB5063877	KB5043050
KB5099538	KB5062557	KB5042215
KB5090354	KB5060531	KB5041578
KB5094123	KB5058392	KB5040430
KB5087061	KB5055519	KB5039879
KB5087538	KB5050533	KB5039217
KB5082123	KB890830	KB5037765
KB5082413	KB5053596	KB5036896
KB5068402	KB5052000	KB5036604
KB5078752	KB5050008	KB5035849
KB5075904	KB5049608	KB5034768
KB5073723	KB5048661	KB5034127
KB5071544	KB5046615	KB5033904
KB5068791	KB5044277	KB5033371
KB5065225	KB5044016	KB5032196
KB5066143	KB5016058	KB4601392
KB5066586	KB5015808	KB4601318
KB5065428	KB5014702	KB4598243

KB5031984	KB5014553	KB4593226
KB5031361	KB5014026	KB4589210
KB5030214	KB5013952	KB4589208
KB5029931	KB5012596	KB4586830
KB5029376	KB5012170	KB4580346
KB5029247	KB5011570	KB4577586
KB5028168	KB5011495	KB4577015
KB5027222	KB5010359	KB4576750
KB5027131	KB5009546	KB4571694
KB5026363	KB5008207	KB4565912
KB5026362	KB5007192	KB4565511
KB5025229	KB5006669	KB4562561
KB5025228	KB5005698	KB4561616
KB5023788	KB5005573	KB4557397
KB5023697	KB5005112	KB4556813
KB5022838	KB5005043	KB4550994
KB5022511	KB5004238	KB4550929
KB5022289	KB5003638	KB4540723
KB5021235	KB5003197	KB4540670
KB5021126	KB5001402	KB4537764
KB5019964	KB5001347	KB4535680
KB5018411	KB5001078	KB4534271
KB5017396	KB5000803	KB4530689
KB5017305	KB5000685	KB4525236
KB5017095	KB4505225	KB4485447
KB5016884	KB4498947	KB4346087
KB5016622	KB4494440	KB4132216

KB4521858	KB4494175	KB4091664
KB4520724	KB4494174	KB2565063
KB4519998	KB4493470	KB2467173
KB4516044	KB4487026	KB4489882
KB4512574	KB4512517	KB4509091
KB4507460		