



Date: August 07, 2026

Author: Software Engineering

System: Cenova

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release July 2026

Purpose:

To release the list of validated Microsoft critical patches on the Cenova System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for July 2026:

KB5100989

KB5099538

KB5090354

List of all patches/KB's validated on the Cenova System:

KB5100989	KB5055519	KB5034127
KB5099538	KB5050533	KB5033904
KB5090354	KB890830	KB5033371
KB5094123	KB5053596	KB5032196
KB5087061	KB5052000	KB5031984
KB5087538	KB5050008	KB5031361
KB5082123	KB5049608	KB5030214
KB5082413	KB5048661	KB5029931
KB5068402	KB5046615	KB5029376
KB5078752	KB5044277	KB5029247
KB5075904	KB5044016	KB5028168
KB5073723	KB5043050	KB5027222
KB5071544	KB5042215	KB5027131
KB5065225	KB5041578	KB5027124
KB5068791	KB5040430	KB5026362
KB5066143	KB5039879	KB5025229
KB5066586	KB5039217	KB5023702
KB5065428	KB5037765	KB5022840
KB5063877	KB5036896	KB5022511
KB5062557	KB5036604	KB5022286
KB5060531	KB5035849	KB5021237
KB5058392	KB5034768	KB5021126

KB5020866	KB5009557	KB5001342
KB5020627	KB5008873	KB5000859
KB5019966	KB5008218	KB5000822
KB5018419	KB5007206	KB5000685
KB5017315	KB5006672	KB4601393
KB5016884	KB5005568	KB4601345
KB5016623	KB5005112	KB4601060
KB5015811	KB5005030	KB4589210
KB5014692	KB5004244	KB4589208
KB5014553	KB5003711	KB4580325
KB5013941	KB5003646	KB4577586
KB5013641	KB5003638	KB4557397
KB5012647	KB5003243	KB4539571
KB5012170	KB5003171	KB4535680
KB5012128	KB5001407	KB4494174
KB5011503	KB5001404	KB4465065
KB5010351	KB5001402	