



Date: July 06, 2026

Author: Software Engineering

System: Cenova

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release June 2026

Purpose:

To release the list of validated Microsoft critical patches on the Cenova System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for June 2026:

KB5094123

List of all patches/KB's validated on the Cenova System:

KB5094123	KB5049608	KB5029247
KB5087061	KB5048661	KB5028168
KB5087538	KB5046615	KB5027222
KB5082123	KB5044277	KB5027131
KB5082413	KB5044016	KB5027124
KB5068402	KB5043050	KB5026362
KB5078752	KB5042215	KB5025229
KB5075904	KB5041578	KB5023702
KB5073723	KB5040430	KB5022840
KB5071544	KB5039879	KB5022511
KB5065225	KB5039217	KB5022286
KB5068791	KB5037765	KB5021237
KB5066143	KB5036896	KB5021126
KB5066586	KB5036604	KB5020866
KB5065428	KB5035849	KB5020627
KB5063877	KB5034768	KB5019966
KB5062557	KB5034127	KB5018419
KB5060531	KB5033904	KB5017315
KB5058392	KB5033371	KB5016884
KB5055519	KB5032196	KB5016623
KB5050533	KB5031984	KB5015811
KB890830	KB5031361	KB5014692
KB5053596	KB5030214	KB5014553
KB5052000	KB5029931	KB5013941
KB5050008	KB5029376	KB5013641

KB5012647	KB5004244	KB4601393
KB5012170	KB5003711	KB4601345
KB5012128	KB5003646	KB4601060
KB5011503	KB5003638	KB4589210
KB5010351	KB5003243	KB4589208
KB5009557	KB5003171	KB4580325
KB5008873	KB5001407	KB4577586
KB5008218	KB5001404	KB4557397
KB5007206	KB5001402	KB4539571
KB5006672	KB5001342	KB4535680
KB5005568	KB5000859	KB4494174
KB5005112	KB5000822	KB4465065
KB5005030	KB5000685	