



Date: September 9, 2026

Author: Software Engineering

System: Faxitron Path Vision

Subsystem: Windows 10 LTSC 2021 & Software Security Patches

Subject: Windows 10 LTSC 2021 & Software Security Patch release August 2026

Purpose:

To release the list of validated Microsoft critical patches on the Faxitron Path Vision System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated:

Date	List of KB's		
July 2026	KB5101006	KB5099539	KB5012170
August 2026	KB5120701	KB5120249	

List of all patches/KB's validated on the Faxitron Path Vision System before July 2026:

KB5101006	KB5094127	KB5012170
KB5099539	KB5087064	KB5082200
KB5012170	KB5087544	KB5082426
KB5066135	KB5078885	