



Date: August 07, 2026

Author: Software Engineering

System: Faxitron Path Vision

Subsystem: Windows 10 LTSC 2021 & Software Security Patches

Subject: Windows 10 LTSC 2021 & Software Security Patch release July 2026

Purpose:

To release the list of validated Microsoft critical patches on the Faxitron Path Vision System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for July 2026:

KB5101006 KB5012170

KB5099539

List of all patches/KB's validated on the Faxitron Path Vision System:

KB5101006 KB5012170

KB5099539 KB5082200

KB5012170 KB5082426

KB5094127 KB5066135

KB5087064 KB5078885

KB5087544