



Date: September 08, 2026

Author: Software Engineering

System: SELENIA DIMENSIONS / 3DIMENSIONS

Subsystem: Windows 10 LTSC 2016 (v1607) & Software Security Patches

Subject: Windows 10 LTSC 2016 (v1607) & Software Security Patch release
August 2026

Purpose:

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated:

<u>Date</u>	<u>List of KB's</u>		
July 2026	KB5099542	KB5099535	KB5090354
August 2026	KB5120418	KB5120236	

List of all patches/KB's validated on the SELENIA DIMENSIONS / 3DIMENSIONS System before July 2026:

KB5099542	KB5065687	KB5042215
KB5099535	KB5063871	KB5041773
KB5090354	KB5062560	KB5041576
KB5094122	KB5062799	KB5040562
KB5094141	KB5060954	KB5040434
KB5087537	KB5061010	KB5039334
KB5088064	KB5058524	KB5039214
KB5082089	KB5058383	KB5037763
KB5082198	KB5055661	KB5037016
KB5068402	KB5055521	KB5036899
KB5078938	KB5050533	KB5035962
KB5075902	KB890830	KB5035855
KB5075999	KB5054006	KB5034862
KB5073447	KB5053594	KB5034767
KB5073722	KB5052006	KB5034119
KB5071543	KB5050109	KB5033373
KB5065225	KB5049993	KB5032391
KB5068864	KB5048671	KB5032197
KB5070247	KB5046612	KB5031362
KB5066584	KB5044293	KB5030504
KB5066836	KB5043124	KB5030213
KB5065427	KB5043051	KB5029376
KB5029242	KB5011570	KB4577586
KB5028169	KB5011495	KB4577015
KB5027219	KB5010359	KB4576750
KB5026363	KB5009546	KB4571694
KB5025228	KB5008877	KB4565912
KB5023788	KB5008207	KB4565511
KB5023697	KB5007192	KB4562561
KB5022838	KB5006669	KB4561616
KB5022289	KB5005698	KB4557397
KB5021235	KB5005573	KB4556813
KB5021126	KB5005043	KB4550994
KB5019964	KB5004238	KB4550929
KB5018411	KB5003638	KB4540723
KB5017396	KB5003197	KB4540670
KB5017305	KB5001402	KB4537764
KB5017095	KB5001347	KB4535680
KB5016884	KB5001078	KB4534271
KB5016622	KB5000803	KB4530689
KB5016058	KB5000685	KB4525236
KB5015808	KB4601392	KB4521858



KB5014702	KB4601318	KB4520724
KB5014553	KB4601051	KB4519998
KB5014026	KB4598243	KB4516044
KB5013952	KB4593226	KB4512574
KB5013625	KB4589210	KB4512517
KB5012596	KB4586830	KB4509091
KB5012170	KB4580346	KB4507460
KB5012118	KB4578969	KB4505225
KB4503537	KB4493470	KB4132216
KB4503267	KB4493448	KB4091664
KB4498947	KB4493435	KB2565063
KB4494440	KB4489882	KB2467173
KB4494175	KB4485447	KB4493472
KB4346087		