



Date: September 08, 2026

Author: Software Engineering

System: SELENIA DIMENSIONS / 3DIMENSIONS

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release August 2026

Purpose:

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated:

<u>Date</u>	<u>List of KB's</u>		
July 2026	KB5100989	KB5090354	KB5099538
August 2026	KB5120698	KB5120238	

List of all patches/KB's validated on the SELENIA DIMENSIONS / 3DIMENSIONS System before July 2026:

KB5100989	KB5099538	KB5090354
KB5094123	KB5050533	KB5034768
KB5087061	KB5049615	KB5034127
KB5087538	KB890830	KB5033904
KB5082123	KB5053596	KB5033371
KB5082413	KB5052000	KB5032196
KB5068402	KB5050008	KB5031990
KB5078752	KB5049608	KB5031361
KB5075904	KB5048661	KB5030214
KB5066143	KB5046615	KB5029925
KB5073723	KB5044277	KB5029376
KB5071544	KB5044016	KB5029247
KB5065225	KB5043050	KB5028168
KB5068791	KB5042215	KB5027222
KB5066137	KB5041578	KB5027124
KB5066586	KB5040430	KB5026362
KB5065428	KB5039879	KB5025229
KB5063877	KB5039217	KB5023702
KB5062557	KB5037765	KB5022840
KB5060531	KB5036896	KB5022504
KB5058392	KB5036604	KB5022286
KB5055519	KB5035849	KB5021237
KB5021126	KB5005112	KB4570333
KB5020874	KB5005030	KB4570332
KB5020615	KB5004244	KB4566424
KB5019966	KB5003711	KB4565349
KB5018419	KB5003646	KB4562562
KB5017315	KB5003243	KB4561608
KB5016884	KB5003171	KB4561600
KB5016623	KB5001407	KB4558998
KB5015811	KB5001404	KB4558997
KB5014692	KB5001342	KB4557397
KB5014553	KB5000859	KB4552924
KB5014164	KB5000822	KB4551853
KB5013941	KB5000685	KB4549949
KB5013641	KB4601393	KB4549947
KB5013626	KB4601345	KB4539571
KB5012647	KB4601060	KB4537759
KB5012170	KB4598480	KB4535680
KB5012119	KB4598230	KB4514366



KB5011503	KB4592440	KB4505225
KB5010351	KB4589208	KB4494174
KB5009557	KB4587735	KB4022619
KB5008873	KB4586793	KB4018073
KB5008218	KB4580325	KB2565063
KB5007206	KB4577668	KB2467173
KB5006672	KB4577667	KB5005568
KB4577586		