



Date: August 07, 2026

Author: Software Engineering

System: SELENIA DIMENSIONS / 3DIMENSIONS

Subsystem: Windows 10 LTSC 2021 & Software Security Patches

Subject: Windows 10 LTSC 2021 & Software Security Patch release July 2026

Purpose:

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for July 2026:

KB5101006	KB5090354
KB5099539	

List of all patches/KB's validated on the SELENIA DIMENSIONS / 3DIMENSIONS System:

KB5101006	KB5087544	KB5075912
KB5099539	KB5082200	KB5073724
KB5090354	KB5082426	KB5071546
KB5094127	KB5068402	
KB5087064	KB5078885	
KB5068781		
KB5065225		
KB5066135		
KB5066791		
KB5065429		
KB5063709		
KB5049613		
KB5012170		
KB5050533		
KB5062554		