



Date: August 07, 2026

Author: Software Engineering

System: SecurView Manager – Windows Server 2016

Subsystem: Windows 10 Server 2016 & Software Security Patches

Subject: Windows 10 Server 2016 & Software Security Patch release July 2026

Purpose:

To release the list of validated Microsoft critical patches on the SecurView Manager – Windows Server 2016 System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for July 2026:

KB5099542

KB5099535

List of all patches/KB's validated on the SecurView Manager – Windows Server 2016 System:

KB5099542

KB5082089

KB5073447

KB5099535

KB5082198

KB5071543

KB5094122

KB5078938

KB5068864

KB5094141

KB5075902

KB5070247

KB5087537

KB5075999

KB5066584

KB5088064

KB5073722

KB5066836

KB5065427	KB5037763	KB4494175
KB5065687	KB5037016	KB4132216
KB5063871	KB5036899	
KB5062560	KB5035962	
KB5062799	KB5035855	
KB5061010	KB5034862	
KB5060954	KB5034767	
KB5058524	KB5034119	
KB5055661	KB5033373	
KB5055521	KB5032391	
KB890830	KB5032197	
KB5054006	KB5031362	
KB5053594	KB5030504	
KB5052006	KB5030213	
KB5050109	KB5029242	
KB5049993	KB5028169	
KB5048671	KB5027219	
KB5046612	KB5026363	
KB5044293	KB5025228	
KB5043124	KB5023788	
KB5043051	KB5023697	
KB5041773	KB5022838	
KB5041576	KB5022289	
KB5040562	KB5017396	
KB5040434	KB5012170	
KB5039334	KB5005698	
KB5039214	KB4589210	