



Date: September 08, 2026

Author: Software Engineering

System: SecurView Manager – Windows Server 2022

Subsystem: Windows 10 Server 2022 & Software Security Patches

Subject: Windows 10 Server 2022 & Software Security Patch release August 2026

Purpose:

To release the list of validated Microsoft critical patches on the SecurView Manager – Windows Server 2022 System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated:

Date	List of KB's		
July 2026	KB5101010	KB5099540	
August 2026	KB5120705	KB5120242	

List of all patches/KB's validated on the SecurView Manager – Windows Server 2022 System before July 2026:

KB5101010	KB5094128	KB5087545
KB5099540	KB5087068	KB5082427
KB5082142	KB5066139	KB5078766