



Date: August 07, 2026

Author: Software Engineering

System: SecurView Manager – Windows Server 2022

Subsystem: Windows 10 Server 2022 & Software Security Patches

Subject: Windows 10 Server 2022 & Software Security Patch release July 2026

Purpose:

To release the list of validated Microsoft critical patches on the SecurView Manager – Windows Server 2022 System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for July 2026:

KB5101010

KB5099540

List of all patches/KB's validated on the SecurView Manager – Windows Server 2022 System:

KB5101010

KB5082427

KB5099540

KB5082142

KB5094128

KB5066139

KB5087068

KB5078766

KB5087545