



Date: August 07, 2026

Author: Software Engineering

System: SecurView

Subsystem: Windows 10 LTSC 2016 (v1607) & Software Security Patches

Subject: Windows 10 LTSC 2016 (v1607) & Software Security Patch release July 2026

Purpose:

To release the list of validated Microsoft critical patches on the SecurView System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for July 2026:

KB5101007	KB5099542	KB5099535
-----------	-----------	-----------

List of all patches/KB's validated on the SecurView System:

	KB5087065	KB5078938
KB5101007	KB5087537	KB5075902
KB5099542	KB5088064	KB5075999
KB5099535	KB5082089	KB5073447
KB5094122	KB5082198	KB5073722
KB5094141	KB5082411	KB5071543

KB5066136	KB5041773	KB5022838
KB5068864	KB5041576	KB5022286
KB5070247	KB5040562	KB5021235
KB5066584	KB5040434	KB5020866
KB5066836	KB5039334	KB5019964
KB5065427	KB5039214	KB5018411
KB5065687	KB5037763	KB5017396
KB5063871	KB5037016	KB5017305
KB5062560	KB5036899	KB5017095
KB5062799	KB5035962	KB5016622
KB5061010	KB5035855	KB5016058
KB5060954	KB5034862	KB5015808
KB5058524	KB5034767	KB5014702
KB5058383	KB5034119	KB5014026
KB5055661	KB5033373	KB5013952
KB5055521	KB5032391	KB5012596
KB890830	KB5032197	KB5012170
KB5054006	KB5031362	KB5011570
KB5053594	KB5030504	KB5011495
KB5052006	KB5030213	KB5010359
KB5050109	KB5029242	KB5009546
KB5049993	KB5028169	KB5008207
KB5048671	KB5027219	KB5007192
KB5046612	KB5026363	KB5006669
KB5044293	KB5025228	KB5005698
KB5043124	KB5023788	KB5005573
KB5043051	KB5023697	KB5005112

KB5005043	KB4561616	KB4346087
KB5004238	KB4561600	KB4132216
KB5003638	KB4557397	KB2565063
KB5003197	KB4556813	KB2467173
KB5001402	KB4550994	
KB5001347	KB4550949	
KB5001078	KB4550929	
KB5000803	KB4540723	
KB4601393	KB4540670	
KB4601392	KB4539571	
KB4601345	KB4537764	
KB4601318	KB4537759	
KB4601060	KB4535680	
KB4598243	KB4534271	
KB4593226	KB4530689	
KB4589210	KB4525236	
KB4589208	KB4521858	
KB4586830	KB4520724	
KB4580346	KB4519998	
KB4580325	KB4516115	
KB4577586	KB4516044	
KB4577015	KB4512574	
KB4576750	KB4512517	
KB4571694	KB4509091	
KB4565912	KB4503308	
KB4565511	KB4494175	
KB4562561	KB4494174	