



Date: September 08, 2026

Author: Software Engineering

System: SecurView

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release August 2026

Purpose:

To release the list of validated Microsoft critical patches on the SecurView System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated:

Date	List of KB's		
July 2026	KB5101008	KB5099538	
August 2026	KB5120703	KB5120238	

List of all patches/KB's validated on the SecurView System before July 2026:

KB5101008	KB5087066	KB5087538
KB5099538	KB5094123	KB5043050
KB5044277	KB5044016	KB5025229
KB5027131	KB5026362	KB5082123
KB5041578	KB5023702	KB5082414
KB5078752	KB5040430	KB5022840
KB5052000	KB5050008	KB5049608
KB5048661	KB5046615	KB5032391
KB5032196	KB5031984	KB5031361

KB5030214	KB5029931	KB5029376
KB5075904	KB5029247	KB5073723
KB5028168	KB5066137	KB5027222
KB5071544	KB5039879	KB5068791
KB5039217	KB5055519	KB890830
KB5053596	KB5035849	KB5034768
KB5034127	KB5034119	KB5033904
KB5033371	KB5066143	KB5066586
KB5063877	KB5065428	KB5037765
KB5062557	KB5060531	KB5058392
KB5036896	KB5036604	KB5022511
KB5022286	KB5021237	KB5021126
KB5020866	KB5020627	KB5019966
KB5018419	KB5017315	KB5016623
KB5015811	KB5014692	KB5013941
KB5013641	KB5012647	KB5012170
KB5012128	KB5011503	KB5010351
KB5009557	KB5008873	KB5008218
KB5007206	KB5005698	KB5005112
KB4589210	KB4589208	KB4577586
KB4535680	KB4494174	KB4132216