



**Date:** July 06, 2026

**Author:** Software Engineering

**System:** SecurView

**Subsystem:** Windows 10 LTSC 2019 (v1809) & Software Security Patches

**Subject:** Windows 10 LTSC 2019 (v1809) & Software Security Patch release June 2026

---

**Purpose:**

To release the list of validated Microsoft critical patches on the SecurView System.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for June 2026:**

KB5094123

**List of all patches/KB's validated on the SecurView System:**

|           |           |           |
|-----------|-----------|-----------|
| KB5094123 | KB5073723 | KB5065428 |
| KB5087066 | KB5066137 | KB5062557 |
| KB5087538 | KB5071544 | KB5060531 |
| KB5082123 | KB5068791 | KB5058392 |
| KB5082414 | KB5066143 | KB5055519 |
| KB5078752 | KB5066586 | KB890830  |
| KB5075904 | KB5063877 | KB5053596 |

|           |           |           |
|-----------|-----------|-----------|
| KB5052000 | KB5031984 | KB5015811 |
| KB5050008 | KB5031361 | KB5014692 |
| KB5049608 | KB5030214 | KB5013941 |
| KB5048661 | KB5029931 | KB5013641 |
| KB5046615 | KB5029376 | KB5012647 |
| KB5044277 | KB5029247 | KB5012170 |
| KB5044016 | KB5028168 | KB5012128 |
| KB5043050 | KB5027222 | KB5011503 |
| KB5041578 | KB5027131 | KB5010351 |
| KB5040430 | KB5026362 | KB5009557 |
| KB5039879 | KB5025229 | KB5008873 |
| KB5039217 | KB5023702 | KB5008218 |
| KB5037765 | KB5022840 | KB5007206 |
| KB5036896 | KB5022511 | KB5005698 |
| KB5036604 | KB5022286 | KB5005112 |
| KB5035849 | KB5021237 | KB4589210 |
| KB5034768 | KB5021126 | KB4589208 |
| KB5034127 | KB5020866 | KB4577586 |
| KB5034119 | KB5020627 | KB4535680 |
| KB5033904 | KB5019966 | KB4494174 |
| KB5033371 | KB5018419 | KB4132216 |
| KB5032391 | KB5017315 |           |
| KB5032196 | KB5016623 |           |