



Date: July 06, 2026

Author: Software Engineering

System: TRIDENT HD

Subsystem: Windows 10 & Software Security Patches

Subject: Windows 10 & Software Security Patch release June 2026

Purpose:

To release the list of validated Microsoft critical patches on the TRIDENT HD System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for June 2026:

KB5094141

KB5094122

List of all patches/KB's validated on the TRIDENT HD System:

KB5094122	KB5088064	KB5082198
KB5094141	KB5082089	KB5068402
KB5087537	KB5082198	KB5078938

KB5075902	KB5042215	KB5021235
KB5075999	KB5041773	KB5021126
KB5073447	KB5041576	KB5021045
KB5073722	KB5040562	KB5019964
KB5071543	KB5040434	KB5018411
KB5065225	KB5039334	KB5017396
KB5068864	KB5039214	KB5017305
KB5070247	KB5037763	KB5017095
KB5066584	KB5037016	KB5016884
KB5066836	KB5036899	KB5016622
KB5065427	KB5035962	KB5016058
KB5065687	KB5035855	KB5015808
KB5063871	KB5034862	KB5014702
KB5062560	KB5034767	KB5014553
KB5062799	KB5034119	KB5014026
KB5060954	KB5033373	KB5013952
KB5061010	KB5032391	KB5012596
KB5055661	KB5032197	KB5012170
KB5055521	KB5031362	KB5011570
KB5050533	KB5030504	KB5011495
KB890830	KB5030213	KB5010359
KB5054006	KB5029376	KB5009546
KB5053594	KB5029242	KB5008207
KB5052006	KB5028169	KB5007192
KB5050109	KB5027219	KB5006669
KB5049993	KB5026363	KB5005698
KB5048671	KB5025228	KB5005573
KB5046612	KB5023788	KB5005043
KB5044293	KB5023697	KB5004238
KB5043124	KB5022838	KB5003638
KB5043051	KB5022289	KB5003197

KB5002464	KB4534271
KB5002406	KB4530689
KB5001402	KB4525236
KB5001347	KB4521858
KB5001078	KB4520724
KB5000803	KB4519998
KB5000685	KB4516044
KB4601392	KB4512574
KB4601318	KB4512517
KB4598243	KB4509091
KB4593226	KB4507460
KB4589210	KB4505225
KB4586830	KB4498947
KB4583465	KB4494440
KB4580346	KB4494175
KB4577586	KB4493470
KB4577015	KB4489882
KB4576750	KB4487026
KB4571694	KB4485447
KB4565912	KB4346087
KB4565511	KB4132216
KB4562561	KB4091664
KB4561616	KB4022619
KB4557397	KB4018073
KB4556813	KB2565063
KB4550994	KB2467173
KB4550929	
KB4540723	
KB4540670	
KB4537764	
KB4535680	