



Date: September 09, 2026

Author: Software Engineering

System: Vision DXA

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release August 2026

Purpose:

To release the list of validated Microsoft critical patches on the Vision DXA System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated:

<u>Date</u>	<u>List of KB's</u>		
July 2026	KB5100989	KB5099538	
August 2026	KB5120698	KB5120238	

List of all patches/KB's validated on the Vision DXA System:

KB5100989	KB5099538	KB5094123
KB5087061	KB5040430	KB5020866
KB5087538	KB5039879	KB5020627
KB5046615	KB5022511	KB5044016
KB5022286	KB5041578	KB5021237
KB5082413	KB5039217	KB5019966
KB5050008	KB5027222	KB5049608

KB5048661	KB5027131	KB5026362
KB5025229	KB5023702	KB5022840
KB5082123	KB5078752	KB4589208
KB5055519	KB890830	KB5032196
KB5053596	KB5031984	KB5052000
KB5031361	KB5030214	KB5029931
KB5029247	KB5028168	KB5075904
KB5073723	KB5071544	KB5068791
KB5066143	KB5037765	KB5066586
KB5060531	KB5058392	KB5035849
KB5034768	KB5034127	KB5033904
KB5065428	KB5063877	KB5036896
KB5062557	KB5036604	KB5018419
KB5017315	KB5016623	KB5015811
KB5013641	KB5012170	KB5005112
KB4577586	KB4535680	