



Date: August 07, 2026

Author: Software Engineering

System: Vision DXA

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release July 2026

Purpose:

To release the list of validated Microsoft critical patches on the Vision DXA System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for July 2026:

KB5100989 KB5099538

List of all patches/KB's validated on the Vision DXA System:

KB5100989	KB5078752	KB5066586
KB5099538	KB4589208	KB5065428
KB5094123	KB5075904	KB5063877
KB5087061	KB5073723	KB5062557
KB5087538	KB5071544	KB5060531
KB5082413	KB5068791	KB5058392
KB5082123	KB5066143	KB5055519

KB890830	KB5027222
KB5053596	KB5027131
KB5052000	KB5026362
KB5050008	KB5025229
KB5049608	KB5023702
KB5048661	KB5022840
KB5046615	KB5022511
KB5044016	KB5022286
KB5041578	KB5021237
KB5040430	KB5020866
KB5039879	KB5020627
KB5039217	KB5019966
KB5037765	KB5018419
KB5036896	KB5017315
KB5036604	KB5016623
KB5035849	KB5015811
KB5034768	KB5013641
KB5034127	KB5012170
KB5033904	KB5005112
KB5033371	KB4577586
KB5032196	KB4535680
KB5031984	
KB5031361	
KB5030214	
KB5029931	
KB5029247	
KB5028168	