

Mycoplasma genitalium Acridinium Ester Labeled Probe Oligo**Analyte Specific Reagent (ASR).**

Analytical and performance characteristics are not established.

Intended Use

Analyte Specific Reagent. Analytical and performance characteristics are not established. This analyte specific reagent is intended to be used as a component of a laboratory developed test (LDT) used exclusively by laboratories. Laboratories can develop their own LDTs using one or more Analyte Specific Reagents and/or general laboratory reagents for *in vitro* diagnostic purposes.

Materials Provided

Mycoplasma genitalium Acridinium Ester Labeled Probe Oligo (Cat. No. PRD-03257)

Component	Quantity	Concentration	Description
Tube	1 x 0.17 mL	0.29 nmol/mL	<i>Acridinium ester labeled probe oligo, complementary to M. genitalium 16S rRNA amplicon.</i>

Warnings and Precautions

- A. For Laboratory Use only. Analytical and performance characteristics are not established.
- B. Use routine laboratory precautions. Wear disposable, powderless gloves, protective eye wear, and laboratory coats when handling specimens and kit reagents. Wash hands thoroughly after handling reagents.
- C. Dispose of all materials that have come in contact with reagents according to local, state, and federal regulations.

Product Specifications**Quality Control**

This analyte specific reagent was manufactured and released in accordance with the Hologic quality control and quality assurance procedures. This reagent has been functionally tested using a system capable of detecting luminescence along with commercially available transcription mediated amplification (TMA) general purpose reagents (GPR). However, it is the responsibility of each laboratory to develop and validate their own tests for clinical diagnostic purposes.

Storage and Stability

Please refer to product labels for information on storage and stability.



Hologic, Inc.
10210 Genetic Center Drive
San Diego, CA 92121 USA

Customer Support: +1 844 Hologic (+1 844 465 6442)
customersupport@hologic.com

Technical Support: +1 888 484 4747
molecularsupport@hologic.com

For more contact information visit www.hologic.com.

Hologic and associated logos are trademarks or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries.

All other trademarks that may appear in this package insert are the property of their respective owners.

This product may be covered by one or more U.S. patents identified at www.hologic.com/patents.

© 2017 Hologic, Inc. All rights reserved.

AW-16922 Rev. 001
2017-07